

Wispa Lite Portable

The Wispa Lite Portable Hoist combines lightweight construction with an integral carry handle enabling easy transportation between different locations. Its simple relocation from room to room negates the need for multiple hoist units.

The Wispa Lite can be used with a full ceiling track installation or free standing gantry system offering controlled transfer from bed to chair; ideal for temporary/short term use as no property adaptation is required.

The compact hoist houses a powerful lift providing controlled lifting for users up to 31 stone (200kg).

Carer Friendly

With its integral ergonomic handle, the Wispa Lite is easy to carry, attach and remove from the ceiling track. Optional attachment handle and lanyards enable comfortable attachment to high ceilings without over reaching or stretching.



User Comfort

Soft start/stop lifting for added user comfort, safety and control.

Enhanced sling attachment design offers wider sling positioning to improve user comfort, whilst providing maximum head clearance.

A room to room transfer system allows the hoist and user to move freely between rooms without the need for repositioning.



The Wispa Lite can be operated using the clearly identifiable buttons on the ergonomically designed handset or by the hoist control panel, thereby reducing risk of not lifting/lowering clients if the handset is lost or damaged.

Visual and audible warnings indicate faults and battery status giving a clear indication of when the hoist should be recharged. The Wispa Lite can simply be plugged into the nearest mains socket and can be fully charged within 2 hours.



Safety

An upper limit switch prevents the lifting tape from over winding into the case and up to the ceiling track.


The slack tape device prevents the lifting tape from winding out of the cassette when an obstruction is encountered during lowering, preventing jamming of the tape.

In the event of hoist failure, the Wispa Lite is fitted with an emergency lowering switch providing safe, controlled descent. An emergency stop with reset facility provides safe operation, should an emergency arise. The hoist is easily reset via the handset or controls on the hoist.




Wispa Lite Portable

Hoist Features



Technical Information



Chiltern Invadex (UK) Ltd
Unit 6C Thorpe Drive
Banbury
Oxfordshire
OX16 4U
Tel: 01869 365500

E-mail: sales@chilterninvadex.co.uk
Web: www.chilterninvadex.co.uk

Overseas enquiries
Tel: 0044 1869 365500
Fax: 0044 1869 365588

Other Product Ranges from Chiltern Invadex

Shower Commode Chairs



TransAqua



AquaMaster A06



AquaMaster Closomat
Range

Tilt in Space Bathing Range



Aquamaster Tilt in Space Chair



Washington Tilting Shower Cradle

Shower Trolleys, Cradles and Stretchers



T1 shower trolley, now available
in 20 colour options



Washington Cradle



T12 fold up stretcher



T3 Paediatric trolley - now
available in 20 colour options

Shower Enclosures, Wet Floors and Cubicles - full installation service available



Shower WC cubicle

EasiDec wet floor

Sandwell level access tray